


NB EOR IN UNCONVENTIONALS

Chip Design Ideas & InspIOR® Configuration



Our Customers (Dec 2024)

World leading energy companies & research organisations



aramco



INPEX



kemira



TENEX



Nanobubble Enhanced Oil Recovery (EOR)

Nanobubbles (50–200 nm) are stable Gas-Filled Particles in Water

- They enhance oil recovery through:

Wettability alteration, interfacial tension (IFT) reduction, penetration into tight pores, capillary pressure modulation

- Why use them?

They're stable, small enough for nanopores, and can carry active gases.

- Potential for “Unconventionals”!

InspIOR & chip design to support testing up to 10.000 psi and 130 °C

Nanobubbles are prepared using a proprietary process available only at UTA facilities.

InspIOR® - Turnkey Microfluidic Technology Platform

Developed for Oil & gas, carbon management, hydrogen, geothermal

- Fast & accurate fluid testing, chemicals screening and flow in porous media analysis
- Specifically designed for high-pressure high-temperature microfluidics applications
- Fully automated system
- Based on 10+ years of microfluidic research & development




What is Microfluidics?

Miniaturized fluid control for faster, smarter lab workflows

- Microfluidics is the science of manipulating tiny amounts of fluids—often nanoliters—through channels thinner than a human hair.

It enables fast, precise, and automated experiments on a compact chip.


CCE / Viscosity / MMP




Dew Point Measurement



MMP – Non Fractured




Phase Behaviour



Chemical Screening



Flow Assurance



Nanobubble EOR - Some Experimental Studies (Only 1 on MF?)

| NO. | REFERENCE & YEAR | INSTITUTION | SYSTEM | METHOD | P RANGE | T RANGE | HIGHLIGHTS |
|-----|------------------------|---------------------------------------|--------------------------------|-----------------------------------|--------------|-----------|--|
| 1 | Shoukry et al. (2025) | University of Wyoming | CO ₂ Nanobubbles | Micromodel (glass, WW/OW) | ~2–6 bar | 25–40 °C | Direct observation of NB transport in high-pressure chips |
| 2 | Cai et al. (2024) | China Univ. of Petroleum (East China) | CO ₂ NB + Pickering | Coreflood (ultra-tight cores) | Up to 10 bar | 25 °C | 66% recovery; pore blocking, wettability shift |
| 3 | Elnaggar et al. (2025) | University of Wyoming | N ₂ Nanobubbles | Core imbibition (Berea/carbonate) | 1–15 bar | 25–60 °C | Enhanced spontaneous imbibition due to NB stability |
| 4 | Li et al. (2024) | Southwest Petroleum Univ. | CO ₂ Nanobubbles | Coreflood (low-perm sandstone) | ~8 bar | 30 °C | Pressure buildup indicates NB mobility control |
| 5 | Wang et al. (2023) | China Univ. of Geosciences | CO ₂ Nanobubbles | Sandpack flooding | ~10 bar | 25–50 °C | Nanobubble-assisted sweep efficiency improved by 20–30% |
| 6 | Zhou et al. (2023) | Univ. of Sci. & Tech. Beijing | CO ₂ Nanobubbles | Coreflood (tight carbonate) | ~5–12 bar | 25–60 °C | NB reduced IFT and improved injectivity |
| 7 | Zhang et al. (2025) | Xi'an Shiyu Univ. / CUP Beijing | CO ₂ Nanobubbles | Coreflood (Berea sandstone) | 5–15 bar | 30–60 °C | Tracked NB diffusion vs CO ₂ gas; improved recovery |
| 8 | Ahmed et al. (2022) | King Fahd Univ. of Petroleum | CO ₂ Nanobubbles | Micromodel (glass, low-perm net) | 3–5 bar | Room temp | Bubble retention and pathway alteration |
| 9 | Wen et al. (2024) | Daqing Oilfield Lab | CO ₂ Nanobubbles | Coreflood (field core, sandstone) | 10 bar | 40–60 °C | NB performance benchmarked vs gas/surfactant |
| 10 | Chen et al. (2023) | China Univ. of Mining & Tech | N ₂ Nanobubbles | Coreflood (tight sandstone) | 5–12 bar | 25–50 °C | High stability N ₂ NB → 10–15% uplift in recovery |


Design 1: Multi-Zone Injectivity Screening

Purpose: Fast & Automated Detection of Pore Entry Threshold of Nanobubbles


- **Chip Design:**
 - Different zones with decreasing channel size (e.g. 1 μm , 0.5 μm , 0.2 μm , 0.1 μm)
 - Each zone has an outlet and microvalve
 - Inlet pressure measured at constant flow rate
- FXL recommends keeping channel width > 5 μm and varying edge depth for better visual access
- **What we learn:**
 - Detect mobility cutoff pore size
 - Identify collapse, retention, or plugging
 - Compare against brine baseline
 - Build injectivity map
 - Vary p,T, salinity, NB concentration, NB size, flow rate (shear rate)

Design 1: Multi-Zone Injectivity Screening

Experimental Setup Sketch



Design 1: Multi-Zone Injectivity Screening (Draft)




FXL recommends keeping channel width $> 5 \mu\text{m}$ and varying edge depth for better visual access

Channel Width

- 0.30 μm
- 0.60 μm
- 0.90 μm
- 1.00 μm
- 2.00 μm
- 3.00 μm

Design 1: Multi-Zone Injectivity Screening (Draft)



Edge Depth

■ 300 nm

FXL recommends keeping channel width $> 5 \mu\text{m}$ and varying edge depth for better visual access

Design 1: Multi-Zone Injectivity Screening

Interpretation Matrix

| ΔP Observation | Possible Interpretation |
|--|--|
| Low ΔP | Nanobubbles freely enter — no blockage or collapse |
| Moderate ΔP increase | Partial retention, bubble-induced flow resistance |
| Sharp ΔP rise at specific zone | Pore entry cutoff reached → plugging or size exclusion |
| Rising ΔP over time (in same zone) | Progressive accumulation → bubble jamming or adhesion |
| Sudden ΔP drop | Bubble collapse or coalescence relieved obstruction |
| High ΔP from Zone 1 onward | Global plugging — bubbles too large or unstable |

NB too small to be visual in microscope.

Design 2: “Knudsen-Analogue” Diffusion of Nanobubbles

Purpose: Fast & Automated Detection of NB Diffusion

- **Bottom channel:** Contains nanobubble-laden brine (source)
- **Top channel:** Contains plain brine (sink)
- **Middle:** Nanochannels
- FXL recommends keeping channel width $> 5 \mu\text{m}$ and varying edge depth for better visual access
- *NB detection via UV-VIS to be tested and calibrated at UTA premises!*
- **What we learn:**
 - Identify the smallest pore size nanobubbles can penetrate.
 - Calculate effective diffusivity from breakthrough time.
 - Detect collapse or retention in nanochannels.
 - Assess injectability in tight formations


Nanobubble – Detection Methods (inline, low volume)


| Method | Principle | Detects | Notes |
|------------------------------------|---|-------------------------|---|
| DLS (Dynamic Light Scattering) | Measures size distribution of particles | Size + presence | Standard for nanobubbles; microfluidic version currently developed by FXL; R&D state |
| UV–Vis Spectroscopy (Transmission) | Measures total extinction (scattering +abs) | Presence (quantitative) | Inline-compatible; expected range 450–600 nm; available at FXL up to ~300 bar; upgrade to 700 bar with nanobubbles to be tested at UTA |
| Conductivity | Non-conductive phase displaces ions | Presence (quantitative) | Δ may be small; affected by ionic strength; needs to be tested; available at FXL; cheap |
| Conv. Microscope | Visual | Size + presence | Challenging to implement; need to test at UTA |

- Other methods not considered due to larger sample volumes required...

Design 2: “Knudsen-Analogue” Diffusion of Nanobubbles

Experimental Setup Sketch – Challenging due to Analytics and dp Requirements



 Micromodel e.g. 200 nm channels

S: Sensor (UV-Vis or DLS) to be tested and calibrated at UTA

Design 2: “Knudsen-Analogue” Diffusion of Nanobubbles

Interpretation Matrix

| Observation | Possible Interpretation |
|----------------------------|-----------------------------------|
| Bubbles appear quickly | Free diffusion through pores |
| Bubbles appear slowly | Hindered diffusion |
| No bubbles appear | Size exclusion or collapse |
| Signal reaches low plateau | Partial transport, some retention |


Design 3: Nanobubble-Based Oil Displacement Micromodel

Visualize and quantify how NB displaces oil in a nanostructure micromodel.

- Inlet / outlet port:
 - Inject / produce oil / brine / NB
- Nanopores:
 - Based on MICP data to resemble real geometrical features of rock
- What we learn:
 - Oil recovery efficiency
 - Displacement mechanism
 - Sweep pattern and front behavior
 - “Residual” oil saturation
 - Effect of pore geometry
- FXL recommends keeping channel width $> 5 \mu\text{m}$ and varying edge depth for better visual access


Design 3: Nanobubble-Based Oil Displacement Micromodel

Experimental Setup Sketch



Micromodel with MICP data-based structure
Homogeneous, vertically or horizontally layered

Design 3: Nanobubble-Based Oil Displacement Micromodel




FXL recommends keeping channel width $> 5 \mu\text{m}$ and varying edge depth for better visual access

Channel Width

- 0.30 μm
- 0.60 μm
- 0.90 μm

Design 3: Nanobubble-Based Oil Displacement Micromodel



FXL recommends keeping channel width $> 5 \mu\text{m}$ and varying edge depth for better visual access.

Continuous path from left to right / (start with conventional chips e.g. sandstone)

May be second step after understanding injectivity. Other data than MICP available? MICP/TS/NMR/FIB-SEM

Design 3: Nanobubble-Based Oil Displacement Micromodel

Interpretation Matrix

| Observation | Possible Interpretation |
|---|---|
| High oil recovery compared to brine | Nanobubbles enhance recovery via wettability or IFT reduction |
| Uniform displacement | Efficient sweep; improved flow |
| Fingering or bypass | Mobility mismatch or channelling; nanobubbles less effective |
| Early breakthrough with low recovery | Poor displacement efficiency; limited bubble function |
| Oil detachment from surfaces | Wettability alteration likely caused by nanobubbles |
| Formation of oil emulsions in effluent | Interfacial disturbance or emulsification during flow |
| Oil remains trapped in corners or dead ends | Insufficient capillary disruption; trapped residual oil |

Way Ahead?

Next Steps and Action Point

- Feedback draft chip designs – anything we missed or should improve?
- What can FXL do to further support process?
- Timeline experimentation at UTA

Adaptive Microfluidics: Systems that learn, optimize & transform labs

fluidXlab GmbH

A Member of the HOT Energy Group

Am Stollen 19B
38640 Goslar, Germany

+49 151 424 407 39

info@fluidxlab.com

www.fluidxlab.com

